

**Patient ID:** 2223314023417  
**Created on:** 11/13/2023  
**Date of birth:** 12/18/2020  
**Receipt of sample:** 11/11/2023  
**Gender:** Male

**Test:** Allergy EUROLINE Inhalation Gulf V2  
**Patient name:** DAYAAN BILAL BILAL  
**Lab number:** 39511  
**Sent from:** PMC ARJAN  
**Age:** 2 years



| Antigen                             | Concentration | Class | 0 | 1 | 2 | 3 | 4 | 5 | 6 |
|-------------------------------------|---------------|-------|---|---|---|---|---|---|---|
| Timothy grass (g6)                  | < 0.35 kU/l   | 0     |   |   |   |   |   |   |   |
| Cultivated rye (g12)                | < 0.35 kU/l   | 0     |   |   |   |   |   |   |   |
| Common ragweed (w1)                 | < 0.35 kU/l   | 0     |   |   |   |   |   |   |   |
| Mugwort (w6)                        | < 0.35 kU/l   | 0     |   |   |   |   |   |   |   |
| Alder (t2)                          | < 0.35 kU/l   | 0     |   |   |   |   |   |   |   |
| Birch (t3)                          | < 0.35 kU/l   | 0     |   |   |   |   |   |   |   |
| Oak (t7)                            | < 0.35 kU/l   | 0     |   |   |   |   |   |   |   |
| Olive tree (t9)                     | < 0.35 kU/l   | 0     |   |   |   |   |   |   |   |
| Dermatophagoides pteronyssinus (d1) | 70.21 kU/l    | 5     |   |   |   |   |   |   |   |
| Dermatophagoides farinae (d2)       | 70.82 kU/l    | 5     |   |   |   |   |   |   |   |
| Cockroach, German (l6)              | < 0.35 kU/l   | 0     |   |   |   |   |   |   |   |
| Cat (e1)                            | < 0.35 kU/l   | 0     |   |   |   |   |   |   |   |
| Dog (e2)                            | < 0.35 kU/l   | 0     |   |   |   |   |   |   |   |
| Horse (e3)                          | < 0.35 kU/l   | 0     |   |   |   |   |   |   |   |
| Penicillium notatum (m1)            | < 0.35 kU/l   | 0     |   |   |   |   |   |   |   |
| Cladosporium herbarium (m2)         | < 0.35 kU/l   | 0     |   |   |   |   |   |   |   |
| Aspergillus fumigatus (m3)          | < 0.35 kU/l   | 0     |   |   |   |   |   |   |   |
| Candida albicans (m5)               | < 0.35 kU/l   | 0     |   |   |   |   |   |   |   |
| Alternaria alternata (m6)           | < 0.35 kU/l   | 0     |   |   |   |   |   |   |   |
| CCD marker (CCD)                    | < 0.35 kU/l   | 0     |   |   |   |   |   |   |   |
| Indicator (Ind)                     | 75.72 kU/l    | 5     |   |   |   |   |   |   |   |
| Label (Et)                          | < 0.35 kU/l   |       |   |   |   |   |   |   |   |

| Concentration (kU/l) | Class | Explanation                                                                                    |
|----------------------|-------|------------------------------------------------------------------------------------------------|
| 0 ≤ sIgE < 0.35      | 0     | No specific antibodies detected.                                                               |
| 0.35 ≤ sIgE < 0.7    | 1     | Very weak antibody detection, often no clinical symptoms with an existing sensitization.       |
| 0.7 ≤ sIgE < 3.5     | 2     | Weak antibody detection, existing sensitization, often with clinical symptoms in upper levels. |
| 3.5 ≤ sIgE < 17.5    | 3     | Definite antibody detection, clinical symptoms are often present.                              |
| 17.5 ≤ sIgE < 50     | 4     | Strong antibody detection, almost always with existing symptoms.                               |
| 50 ≤ sIgE < 100      | 5     | Very high antibody titer.                                                                      |
| sIgE ≥ 100           | 6     | Very high antibody titer.                                                                      |

