

Patient ID: 2224307028617

Created on: 11/4/2024

Date of birth: 1/1/1972

Receipt of sample: 11/2/2024

Gender: Male

Test:

Patient name:

Lab number:

Sent from:

Age:

Allergy EUROLINE Food Gulf\_V3

DHEERAJ MAHENDRU

26316

CITICARE MEDICAL CENTER

52 years



| Antigen              | Concentration | Class | 0 | 1 | 2 | 3 | 4 | 5 | 6 |
|----------------------|---------------|-------|---|---|---|---|---|---|---|
| Egg white (f1)       | < 0.35 kU/l   | 0     |   |   |   |   |   |   |   |
| Egg yolk (f75)       | < 0.35 kU/l   | 0     |   |   |   |   |   |   |   |
| Cow's milk (f2)      | < 0.35 kU/l   | 0     |   |   |   |   |   |   |   |
| Wheat flour (f4)     | < 0.35 kU/l   | 0     |   |   |   |   |   |   |   |
| Soybean (f14)        | < 0.35 kU/l   | 0     |   |   |   |   |   |   |   |
| Baker's yeast (f45)  | < 0.35 kU/l   | 0     |   |   |   |   |   |   |   |
| Nut mix 2 (fs36)     | < 0.35 kU/l   | 0     |   |   |   |   |   |   |   |
| Peanut (f13)         | < 0.35 kU/l   | 0     |   |   |   |   |   |   |   |
| Orange (f33)         | < 0.35 kU/l   | 0     |   |   |   |   |   |   |   |
| Strawberry (f44)     | < 0.35 kU/l   | 0     |   |   |   |   |   |   |   |
| Banana (f92)         | < 0.35 kU/l   | 0     |   |   |   |   |   |   |   |
| Mango (f93)          | < 0.35 kU/l   | 0     |   |   |   |   |   |   |   |
| Tomato (f25)         | < 0.35 kU/l   | 0     |   |   |   |   |   |   |   |
| Carrot (f31)         | < 0.35 kU/l   | 0     |   |   |   |   |   |   |   |
| Onion (f48)          | < 0.35 kU/l   | 0     |   |   |   |   |   |   |   |
| Chicken (f83)        | 0.35 kU/l     | 1     |   |   |   |   |   |   |   |
| Lamb meat (f88)      | < 0.35 kU/l   | 0     |   |   |   |   |   |   |   |
| Codfish (f3)         | < 0.35 kU/l   | 0     |   |   |   |   |   |   |   |
| Shrimp/ Prawn (f24)  | < 0.35 kU/l   | 0     |   |   |   |   |   |   |   |
| CCD marker (CCD)     | < 0.35 kU/l   | 0     |   |   |   |   |   |   |   |
| Indicator band (Ind) | >= 100 kU/l   | 6     |   |   |   |   |   |   |   |
| Label (Et)           | < 0.35 kU/l   |       |   |   |   |   |   |   |   |

| Concentration (kU/l)          | Class | Explanation                                                                                    |
|-------------------------------|-------|------------------------------------------------------------------------------------------------|
| $0 \leq \text{slgE} < 0.35$   | 0     | No specific antibodies detected.                                                               |
| $0.35 \leq \text{slgE} < 0.7$ | 1     | Very weak antibody detection, often no clinical symptoms with an existing sensitization.       |
| $0.7 \leq \text{slgE} < 3.5$  | 2     | Weak antibody detection, existing sensitization, often with clinical symptoms in upper levels. |
| $3.5 \leq \text{slgE} < 17.5$ | 3     | Definite antibody detection, clinical symptoms are often present.                              |
| $17.5 \leq \text{slgE} < 50$  | 4     | Strong antibody detection, almost always with existing symptoms.                               |
| $50 \leq \text{slgE} < 100$   | 5     | Very high antibody titer.                                                                      |
| $\text{slgE} \geq 100$        | 6     | Very high antibody titer.                                                                      |

| Band      | Band Information                  |
|-----------|-----------------------------------|
| Nut mix 2 | Contains Peanut, Hazelnut, Almond |

